**AquaVial™ Plus** is designed for testing well and ground water and consists of two tests: **AquaVial™ Total Bacteria** and **AquaVial™ E. coli**. **AquaVial™ Plus** is a quick and convenient screening test and does not replace mandatory water testing through an accredited laboratory.

**AquaVial™ Total Bacteria** detects bacteria and fungi concentrations as low as 500 CFU/ml (500 bacteria per ml) in 15 minutes at 35-40 °C (95-104 °F), or in 30 minutes at 20-25 °C (68-77 °F).

**AquaVial™ E. coli** detects E.coli and coliform concentrations as low as 1 CFU/ml (one bacteria per ml) of E. coli and coliform bacteria within 24 hours when incubated at 35-40 °C (95-104 F), or 48 hours at 20-25 °C (68-77 °F).

#### KIT CONTENT:

- 1 T-Cup unit
- 1 AquaVial™ Total Bacteria vial with **pink solution**
- 1 AquaVial™ E. coli vial with **yellow dry media**
- 2 AquaVial™ plug caps
- 1 distilled water pouch
- 1 set of Testing Instructions

#### IMPORTANT SAFETY NOTES:

- The solutions in the vials are non-toxic and; therefore pose no significant hazard provided that the directions for use are followed. **DO NOT INGEST.**
- Not recommended for children (small parts)
- Store at temperatures above freezing. For best results, test before the expiration date indicated on the box

#### TIP:

Have a white background (i.e. a piece of paper) and a camera ready before starting the tests. This will help observing the colour change.

#### DISPOSAL:

T-Cup components, are recyclable. Please recycle whenever possible. Vials can be disposed of in the trash.

**Warning: Do not attempt to open the vials after test completion or dispose of their contents into your sink.**

#### IMPORTANT:

**Please read the directions carefully and watch the instructional video, before starting the test.** Follow the steps below in the exact order to obtain an accurate reading of your water quality. Any deviation may cause cross contamination or ruin the sample, possibly resulting in false positive or false negative results

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
##### Step 1: Collect and Concentrate Water Sample

**OPEN T-CUP:** Separate the top and bottom halves of the T-Cup by twisting the top half **counter clockwise**.

**ADD WATER SAMPLE:** Fill the bottom of the T-Cup **halfway** with the sample water.

**CONCENTRATE:** Twist top of the T-Cup **clockwise** into bottom section until it stops, forcing water through the integrated T-Cup filter. Water will pour through and spill over the sides of the plunger. **This is normal.** The T-Cup is intended to sample this exact amount of water.

**DISCARD:** Drain excess water from the top and sides of the T-Cup.

##### Step 2: Add Distilled Water


This step is required to ensure the water sample collected and concentrated in Step 1 is washed of chemicals, like chlorine, used to treat the water, which may affect your test results. The addition of distilled water will not dilute or interfere with the water sample to be tested. The T-Cup system has a built-in filter that retains the bacteria at the bottom, while eliminating chemicals through the filter.

**OPEN T-CUP:** Open the T-Cup by twisting the top half counter-clockwise.

**ADD DISTILLED WATER:** Carefully tear the top section of the distilled water pouch and pour distilled water into the bottom half of the T-Cup.

**CONCENTRATE:** Twist the top of the T-Cup into the bottom clockwise to close the T-Cup.

**DISCARD:** Pour out any excess water from the top of the T-Cup.

##### Step 3: Test Water Sample

**OPEN T-CUP:** Open the T-Cup by twisting the top half counter-clockwise.

**SNAP OPEN VIAL:** Hold AquaVial™ Total Bacteria vial (pink solution) into your hand **firmly** in a vertical position with the tip of the cap just above your index finger. Apply side pressure with your thumb to snap the cap open. **BE CAREFUL NOT TO TILT OR SPILL AQUAVIAL™**

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**FILL AQUAVIAL™:** Pour the unfiltered water from the bottom of the T-Cup into the AquaVial™ Total Bacteria vial. In most cases the liquid will not reach the fill line marked on the vial, which is normal. Sample water should not exceed the fill line; or stop filling water in the vial after this line. Close AquaVial™ Total Bacteria vial using one of the two plug caps provided.


**Tip:** Take a picture of the vial against a white background.

**INCUBATE:** Let the AquaVial™ Total Bacteria vial incubate in an upright position for 30 minutes and check for changes in colour. If you are not sure if the vial changed color, let the vial incubate another 30 minutes and check color again.

Compare the color of the reagent in the test vial against the color diagram below:

35-40°C (95-104°F) – 15 min OR 20-25°C (68-77°F) – 30 min

If the solution changes its colour to clear, **extreme caution should be exercised**, and a sample should be taken to a laboratory for further testing. This test is intended for **rapid screening only**.

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#### Step 1: Test Water Sample


**SNAP OPEN VIAL:** Hold AquaVial™ E. coli vial (yellow dry media) into your hand **firmly** in a vertical position with the tip of the cap just above your index finger. Apply side pressure with your thumb to snap the cap open.

**FILL AQUAVIAL™:** Pour sample water directly from the source into the AquaVial™ E. coli vial. Sample water should not exceed the fill line; or stop filling water in the vial after this line. Close AquaVial™ E. coli vial using one of the two plug caps provided. **DO NOT LEAVE THE VIAL OPEN.**

**SHAKE:** Shake the vial well to dissolve the dry media deposited on the side wall of the vial.


**Tip:** Take a picture of the vial against a white background.

**INCUBATE:** Let the AquaVial™ E. coli vial incubate in an upright position for 24 hours at 35 - 40°C (95 - 104°F) or 48 hours at 20-25°C (68 - 77°F) and check for changes in colour. Compare the color of the reagent in the test vial against the colour diagram below:

Compare the color of the reagent in the test vial against the color chart below:

35-40°C (95-104°F) – 24 hrs OR 20-25°C (68-77°F) – 48 hrs

If the solution changes its colour to orange or purple, extreme caution should be exercised, and you should not come into contact, ingest or shower with this water, until further testing has been conducted and additional analysis performed. This test is intended for rapid screening only.

#### RESULTS INTERPRETATION:

Untreated well, ground of surface water normally contains bacteria, fungi and viruses.

High levels of total bacteria and fungi (over 500 CFU/ml) are an indicator of untreated water. Very high bacteria levels (over 2000 CFU/ml) indicate biofilm presence. Biofilm is a slimy, sticky film, produced by bacteria as a way to store nutrients and defend against chemicals such as Chlorine. The presence of biofilm in warm waters (in cooling towers, hot tubs, showers) has been associated with the presence of disease – causing bacteria such as Legionella, Pseudomonas, and Klebsiella - bacteria known to cause Pneumonia.

E coli and coliform presence is an indicator of water being contaminated with wild or farmed animal feces, and can cause pose serious health risks.

World Health Organization recommends that drinking and recreational water have Less than 500 CFU/ml of total bacteria and fungi and less than 1 CFU/ml of E coli and coliform bacteria.

**A test result of less than 500 CFU/ml for total bacteria and fungi and no E coli and coliform is a very good indicator that the water is safe and properly treated. Any results over these limits should be confirmed by a certified lab as soon as possible, and the water should not be used until the results are received from the laboratory confirmatory tests. Alternatively, an effective water treatment should immediately be implemented.**

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